UNIVERS MAV

Recirculating Fume Hood Adsorption on Activated Carbon Filter



The Recirculating Fume Hood, type UNIVERS MAV, is intended to protect the user and its environment close to volatile, harmful or malodorous products thanks to a system of vertical aspiration on activated carbon.

Double Protection:

Operator Environment

Complies with standard:

- Self-assembly structure
- Manual regulation
- Glass volume
- tilting frontal window



	Univers MAV 6		Univers MAV 9	
	internal	Overall	internal	Overall
Width (mm)	600	620	905	925
Height (mm)	905	925	905	925
Depth (mm)	592	600	592	600
Weight of liter of CA	25		30	