





Protection : Operator

✓ No maintenance✓ Custom made

TPL version with ventilated cupboard

The T is designed to protect the user and its surrounding environment by means of an exhaust system. Made of non-particle-generating materials, they can be installed directly in clean rooms.

The optimization of the dimensions makes it possible to integrate in any type of laboratory, even low ceiling.

	T 81		T 82		T 83		T 84	
	internal	Overall	internal	Overall	internal	Overall	internal	Overall
Width	784	800	784	800	1204	1220	1484	1500
Height	267	275	internal: 642 OVERALL: 650					
Depth	internal: 340 OVERALL: 590							
Diameter	125		160		200		200	
Flow rate (m^3/h)	300		700		1000		1200	