

PPS

Biosafety Cabinet Type II



Triple protection :
Operator
Product
Environment

Complies with standard:
NF EN 12469

- ✔ Electric sash window
- ✔ Bag in Bag out System
- ✔ Maintenance from the front
- ✔ 3 HEPA filters

★ FLAGSHIP PRODUCT

Microbiological safety stations type II are designed to protect the operator, product and the environment from any airborne contamination.

The PPS is suitable for handling pathogenic germs of class 1, 2 and 3.
The PPS can be fitted with a marble for accurate weighing at (0,01 mg)



	PPS 9	PPS 12	PPS 16	PPS 18
Useful Width	921	1226	1531	1836
Total width	971	1276	1581	1886
Height	internal : 650		overall : 2200	
Depth	Useful: 616		Overall: 786	
Weight (kg)	270	330	420	490
Power (KW)	1,4	1,4	1,5	1,5