

Air Dépoussiéré Stérile Laminaire

# TECHNICAL FEATURES

## MINISPACE

Protection of the Manipulation

Laminar Flow cabinet ISO 5

Comply with standard: ISO 14644









catalogues

#### Solutions de l'Ultra Propreté

34, rue Maurice de Broglie - 93600 - Autnoy Sous Pois - France

tél:+33 1 48 17 85 10 / Fex:+33 1 48 17 85 17

Web: www.adslaminaire.cem / E-mail: info@ledslaminaire.com

#### **MINISPACE**



PROTECTION OF THE MANIPULATION

• Automatic ECM motor fan

• HEPA filter H14 99.999%

Led lighting

Transparent work area

Stainless steel work surface


• Prefilter

Control panel with velocity dispaly

Flip window for total opening

Ergonomic inclined front window





ISO Class 5 laminar flow units are intended to protect . manipulation through a vertical ultra-clean airflow

The MINISPACE's blowing and filtration system guarantees work in a sterile environment.

Fully automatic, the MINISPACE has LED lighting and is suitable for all laboratory configurations.

	MINISPACE	
	internal	overall
width(mm)	780	800
height (mm)	internal: 650 – Hors Tout: 935 on stand 1766	
depth (mm)	Utile: 500 - overall: 510	
weight (Kg)	50	



the following documents are delivered:

- certificate of conformity
- warranty certificate
- user's manuel





### Solutions de l'Ultra Propreté

34, rue Maurice de Broglie - 93600 - Autna / Sins Bois - France tél : + 33 1 48 17 85 10 / Fex : +33 1 48 17 85 17

Web: www.adslaminaire.com / E-meil: info@delslaminaire.com