IV GEL SPACE

Vertical laminar flow - Class ISO 5





Complies with standard:

- **⊘** ECM regulation
- **✓** Full front opening
- **⊘** Glass volume



The enclosures Type IV GEL consist of a blowing and filtration unit, guaranteeing work in a sterile environment.

The cabinets class ISO 5 are designed to protect the handling thanks to a vertical air flow.

The filtration system based on a self-regulating unit type FANJET, guarantees a work in sterile and clean environment of ISO 5.

	IV GEL SPACE 9		IV GEL SPACE 12		IV GEL SPACE 16		IV GEL SPACE 18	
	internal	Overall	internal	Overall	internal	Overall	internal	Overall
Width (mm)	885	930	1190	1223	1556	1573	1792	1825
Height (mm)	internal: 650 OVERALL: 1015							
Depth (mm)	internal : 590 OVERALL : 600							
Height base(mm)	1870							
Weight (Kg)	80		10	00	130		150	